

# DATASHEET

## CamCover Sterile Video Camera Light Handle Cover

Code Number/s: TC65470

Our single-use, sterile light handle cover is designed for use in conjunction with light handles with built-in video cameras:

- Foam ring which holds the cover in place
- Clear lens to ensure no video interference
- Latex free
- Peel-off cover to ensure no scratching of lens

### Product Image



### Product Information

<b>Product Description</b>	Sterile Video Light Handle Cover
<b>Proposed Use</b>	To help protect the sterile field by providing a sterilised cover for the video light handle.
<b>Packing Quantity</b>	120 per box Individually Sterile Packed
<b>Sterilization</b>	Ethylene Oxide
<b>Shelf Life</b>	5 Year shelf life
<b>Compliance</b>	Medical Devices Directive 93/42/EEC
<b>Class</b>	1
<b>Manufacturer</b>	Shanghai Everpure Medical

# Toffeln<sup>®</sup>

## Packaging

### Labelling of product

Each unit

- Article number,
- Lot number,
- Manufacturer
- Production date,
- Expiry date,
- Single use
- CE

### Storage Instructions

#### Handling and Storage:

Product should be stored in its original closed carton, in a dry place.

Do not expose to humidity and direct heat or sunlight.