

TC65SCWM ViewMax Electric Scope Warmer

The ViewMax Electric Scope Warmer TC65SCWM is a versatile, cost-effective and reliable solution that maintains clear visibility during surgery by preventing and eliminating condensation on the laparoscope lens.

Rechargeable Battery Pack TC65SCWB sold separately, but included in first order.



Features

- » **Pre-Sterilised Warming Probe:** Heats the laparoscope's distal lens to body temperature.
- » **Adjustable Temperature Control:** Raise, decrease and holds the temperature at 4 different setting between 35°C - 65 °C.
- » **Rechargeable Battery and Adapter:** Rechargeable battery and adapter included on first order and available to purchase separately thereafter.
- » **Compliance:** Product is CE marked and compliant to Medical Device Directive 93/42.

Benefits

- » **Consistently Clear Visibility:** Eliminates fogged lenses during lengthy procedures. The scope is heated to body temperature, eliminating the need for gels or wipes and ensuring optimal clarity.
- » **Enhanced Surgical Efficiency:** Customisable temperature settings allow you to fine-tune the device to your needs, reducing interruptions and maintaining workflow.
- » **Ready to Use:** Each scope warmer comes pre-sterilised for immediate use, and the rechargeable battery can be recharged after each procedure.

ViewMax Scope Warmer



ViewMax Rechargeable Battery & Adapter



Product Specification

ViewMax Scope Warmer

- » **Scope Warmer Dimensions:** 200mm long, designed for up to 12mm in diameter laparoscopes.
- » **Sterile:** Pre-sealed and sterilised, and ready for immediate use.

Rechargeable Battery & Adapter

- » **Temperature Control:** Adjusts temperature via the rechargeable battery pack with a simple click system:
 - » **1 Click** Red Light ● : 65°C
 - » **2 Clicks** Orange Light ● : 55°C
 - » **3 Clicks** Green ● : 45°C
 - » **4 Clicks** Flashing Green ● : 65°C
 - » **5 Clicks** Off ●
- » **Battery Type:** Smart Lithium-ion
- » **Voltage:** 7.4V
- » **Capacity:** 4400mAh
- » **Input:** DC 8.4V, 1000mA
- » **Output Power:** Maximum 15W
- » **Adapter:** For battery pack recharging.



TC65SCWM



TC65SCWB



Packaging & Storage

Dimensions (LxWxD) & Weight per pack:

Scope Warmer:	300 x 100 x 20 (mm) 50g
Rechargeable Battery Pack	100 x 80 x 20 (mm) 220g
Adapter:	100 x 70 x 55 (mm) 90g

Scope Warmer Storage: Store in a dry environment away from direct sunlight. If stored correctly, the Electric Scope Warmer has a 3-year shelf life from the date of manufacture.

Ordering Information

Product	Product Code	Pack Weight	Qty
Scope Warmer	TC65SCWM	1.3kg	25
Battery Pack & Adapter*	TC65SCWB	310g	1
Adapter	TC65SCWC-T	90g	1

*Included free of charge on first order.