



## TG015

### Disposable Short Sleeve Warming Jacket

75% Recycled material


The Disposable warming jacket is made from 75% recycled spunbonded fabric which is soft and breathable. The perfect solution to cold operating theatres. The Toffeln Warming Jacket ensures that theatre staff are kept warm and comfortable throughout the day. Being single-use this helps manage effective infection control and reduces the risk of microorganism contamination.

#### Features

- › Extra layer of comfort
- › Compatible with infection control measures
- › Manufactured from recycled low lint, comfort fabric
- › Ultrasonic welded seams
- › Knitted cuffs
- › Popper stud front



#### Materials



##### EcoTex

Manufactured from a mix of 75% Recycled SS spunbonded and SMS spunbonded melt-blown fabric. It is soft and comfortable to wear and is embossed with a texture to provide no-see through modesty properties. It is a favourite for patient gowns, scrubs and warming jackets.

#### Quality Assurance


Our Disposable Warming Jacket is manufactured in accordance with the EN ISO 13485 standard. This ensures that we deliver products and services that consistently meet customer requirements, guaranteeing full customer satisfaction.



### Sizing (CM)

Medium TG015M  
 Large TG015L  
 Extra Large TG015XL  
 Extra Extra Large TG015XXL

SIZE	M	L	XL	XXL
<b>A</b> Chest	67	70	72	75
<b>B</b> Length	77	80	84	87
<b>C</b> Hem Width	65	68	70	73
<b>D</b> Sleeve (Full)	74	76	78	80
<b>E</b> Neck Opening	17	19	20	21
<b>F</b> Cuff Length	8	8	8	8



### Packaging & Storage Details

Store at ambient temperature in dry conditions.  
 The product has a shelf life of 5 years.

**Box Quantity** 50 pcs

**Dimensions** 44 x 40 x 60 cm